



GIRAFFE

BOTTLE

CARE & USE GUIDE



Introduction

The Giraffe Bottle™ Hands Free Drinking Bottle allows users with various abilities to stay hydrated. The system is designed to be flexible and easy to use. The Giraffe Bottle™ is a great solution for desks, wheelchairs, beds, tables, and more.

The Giraffe Bottle™ is compatible with our aluminum bottle cage, weighted desktop base, and a range of mounting solutions offered by MOGO Mounts.

Contents

BOTTLE

750ml/25oz Tritan Bottle. FDA compliant BPA free plastic; it is not manufactured with bisphenol A (BPA) or other bisphenol compounds, such as bisphenol S (BPS)

LID WITH SEAL & CAP

FDA compliant polypropylene plastics with silicone seal

NECK

FDA compliant material acetal copolymer

DRINKING TUBE


BPA and phthalates free
FDA compliant

CHECK VALVE

FDA compliant polypropylene with silicone seal

CLEANING BRUSH

Stainless steel with FDA compliant nylon brush



Maximum Operating Temperature: 185°F / 85°C

Assembly Guide

This drinking solution is designed to be quickly disassembled for easy cleaning, travel, and routine part replacement.

It is recommended that you disassemble and clean prior to first use.


INITIAL ASSEMBLY:

- 1.) Screw the neck assembly onto the lid base.
- 2.) Install the lid at the top of the bottle.

FULL DISASSEMBLY:

- 1.) Lid unscrews from the bottle.
- 2.) The flexible neck assembly unscrews from the spout of the lid.
- 3.) Uninstall the check valve from the drinking tube.
- 4.) Remove the drinking tube - straighten out the neck and pull tube out at one end.
- 5.) The lid seal can be removed and cleaned.

For assembly, reverse the disassembly procedure.



TIPS

To help reinsert the drinking tube into the neck, roll the neck assembly on a table to straighten the beads. Inserting the tube by feeding it from the bottom will make reassembly easier.

Neck length can be adjusted if necessary. Every plastic bead can be removed if enough force is applied. An assembly tool, available separately, can be used to assist in reassembling the neck. However, assembly is possible by hand.

Cleaning

Thoroughly wash after each use.

For best results, disassemble and place the bottle, bottle lid, and lid seal in a dishwasher. Hand wash the drinking tube and check valve in a sink/container of warm soapy water. The neck assembly must be hand-washed independent of other components.


DISHWASHER

While all components are top-shelf dishwasher safe, it is recommended you only use the dishwasher for the bottle, bottle lid, and lid seal. Some detergents contain ingredients that may discolor or shorten the life of components. It is recommended to limit the use of detergents with bleach.

HAND-WASH



It is recommended that you hand-wash the drinking tube and check valve in a sink or container. For best results, submerge the drinking tube in soap and hot water and use the supplied cleaning brush. Rinse all components and leave to dry before use.



When cleaning the complete neck assembly, **do not** submerge in water. To clean, gently wipe the spout tip and threaded base with a damp cloth. Use the supplied cleaning brush to gently scrub the inlets of both ends and wipe dry.

Additional cleaning of the neck assembly may be required periodically, depending on the user. However, due to difficulty in reassembly, we do not recommend complete disassembly for reasons other than thorough, deep cleanings.



Replacement Schedule

Following regular, daily cleaning, it is still recommended that integral parts of the Giraffe Bottle™ Hands Free Drinking System be replaced periodically to maintain optimal bottle hygiene and usability.

REPLACEMENT PART SCHEDULE

The components that require routine replacement are detailed below. You should replace each part in accordance with the replacement schedule to achieve the best experience with the Giraffe Bottle™ Hands Free Drinking System. The schedule reflects daily use accompanied by routine cleanings.

DRINKING TUBE

MONTHLY

CHECK VALVE

MONTHLY

CLEANING BRUSH

EVERY 3 MONTHS

NECK ASSEMBLY

EVERY 8 MONTHS

Your neck assembly replacement schedule will vary depending on the length configuration of your neck and your ability to keep the whole assembly clean. We recommend replacing the neck, on average, every 8 months. If your bottle has a longer neck, it might require replacement more frequently for bottle hygiene and function.

Alternate Uses

The Giraffe Bottle™ Hands Free Drinking System was designed for those with assistive needs.

However, your Giraffe Bottle can be used as a direct-drink hydration bottle. No matter your needs, Giraffe Bottle has “exactly what Gir-after”.

Nondescript orange fluid
not included. →

More info at
GiraffeBottle.com





DRINKING TUBE & CHECK VALVE

Drinking tube with check valve for monthly replacement.

#504082 for 15" neck

#504084 for 25" neck

#504085 for 30" neck

#504086 for 36" neck

These parts should be replaced on a monthly basis; replacement schedule will vary depending on use.



DRINKING TUBE

Replacement drinking tube.

#504512 24" tube for 15" neck

#504514 36" tube for 25" neck

#504515 41" tube for 30" neck

#504516 45" tube for 36" neck



CLEANING BRUSH

25" stainless steel stem with nylon brush.

#504711



CHECK VALVE

Replacement check valve.

#504611



NECK ASSEMBLY

Flexible, plastic neck assembly. Does not include drinking tube.

#504012 15" neck assembly

#504014 25" neck assembly

#504015 30" neck assembly

#504016 36" neck assembly



SUCTION BULB

Medical-grade suction bulb allows the fluid to be brought to the end of the drinking tube.

#504711



BOTTLE

750ml/25oz Tritan Bottle with lid, cap, and seal. This part can be used as a replacement or as an additional bottle.
#504111




TALL BOTTLE JACKET

Put a jacket on your giraffe to keep it warm or cool.
#504053



BOTTLE HOLDER

Aluminum bottle holder, clamp rotates horizontally & vertically.
#500311



DESKTOP BOTTLE HOLDER

Weighted base for holding bottles on a table, desk, or other flat surface. Includes a Bottle Jacket.
#500322



Recommendations

- Clean after every use.
- Keep the bottle lower than the tip of the straw. If the bottle is higher than the tip of the straw it can create a siphon effect. Once liquid begins to flow down, gravity causes it to continue flowing.
- Do not microwave.
- For best results, keep the neck as short as your application will allow.
- Not recommended for hot beverages such as coffee or hot chocolate; the bottle could become too hot to touch.
- Water and thin liquids work best.
- The check valve does not work with carbonated beverages, drinks with pulp, or other thick liquids. Simply remove the check valve when enjoying these types of beverages.
- With the attached cap, the adjustable neck and drinking tube can easily be removed and the bottle can be put into the refrigerator or exchanged with another bottle.
- If the drinking tube is too flexible at the end, try shortening the tube that extends past the neck. However, it is recommended that you keep the drinking tube extended beyond the neck; refrain from sipping on the end of the neck assembly.
- Do not freeze check valve in liquid.
- Do not turn the bottle over with the neck attached, it will leak.

How to Order

You can visit us online and browse our full catalog of Giraffe Bottle products. While all items can be ordered online, you can also submit your order by phone or email.

WWW.GIRAFFEBOTTLE.COM
SALES@GIRAFFEBOTTLE.COM



GIRAFFE
B O T T L E

A DIVISION OF:

 **Oak Hill**
BRANDS CORP

Product designed in the USA.

Manufactured with domestic and foreign components.

Information and specifications subject to change at any time without notice.

Giraffe Bottle is a registered trademark of Oak Hill Brands Corp.